


Anti-Atpb Mouse Monoclonal Antibody For Wb, Ihc-P, Elisa - P06576

Item Number: CM0006787



Introduction

High-specificity mouse monoclonal antibody against human ATP5B. Validated for WB, IF/ICC, IHC-P, and ELISA; cross-reacts with mouse and rat. Suitable for ATP synthase beta subunit and mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	ATPB Mouse mAb
Remarks/Alias	ATP5B; ATPMB; ATPSB; HUMOP2; HEL-S-271
Species	Human
GeneID (Human)	506
GeneID	506
Immunogen	Synthetic peptide - A synthesized peptide derived from human ATPB
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06576
Protein Weight	57kDa
Shipping	Ice bag