


Anti-Atpb Rabbit Polyclonal Antibody - P06576

Item Number: CM0030274



Introduction

Rabbit polyclonal antibody targeting mitochondrial ATP synthase subunit beta (ATPB/ATP5B). Validated for WB, IHC-P, IF/ICC, and ELISA with reactivity in human, mouse, and rat. Ideal for mitochondrial energy metabolism research.

[Learn More](#)

Parameter	Value
Product Name	ATPB Rabbit pAb
Remarks/Alias	ATP5B; ATPMB; ATP5B; HUMOP2; HEL-S-271; ATPB
Species	Human
GeneID (Human)	506
GeneID	506
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 230-529 of human ATPB (NP_001677.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06576
Protein Weight	57kDa
Shipping	Ice bag