


Anti-Atpif1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Uii2

Item Number: CM0005362



Introduction

Rabbit monoclonal anti-ATPIF1 antibody validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat ATPIF1. Ideal for mitochondrial function, ATPase inhibition, and heme synthesis studies.

[Learn More](#)

Parameter	Value
Product Name	ATPIF1 Rabbit mAb
Remarks / Alias	IP; ATP1; ATPIP; ATPIF1
Species	Human
Gene ID (Human)	93974
Gene ID	93974
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-106 of human ATPIF1 (NP_057395.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Uii2
Protein Weight	12kDa
Shipping	Ice bag