


Anti-Atpif1 Rabbit Polyclonal Antibody (Ko Validated) For Wb And Ip - Q9Uii2

Item Number: CM0030942



Introduction

Rabbit polyclonal antibody against human Atpif1, validated for knockout (KO) specificity. Suitable for WB, IP, IF/ICC, ELISA. Cross-reacts with mouse. Ideal for mitochondrial ATPase inhibition and heme synthesis research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Atpif1 Rabbit pAb
Remarks/Alias	IP; Atpi; Atpip; Atpif1; F1
Species	Human
Gene ID (Human)	93974
Gene ID	93974
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-106 of human Atpif1 (NP_057395.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Uii2
Protein Weight	12kDa
Shipping	Ice bag