


Anti-Atr Polyclonal Antibody For Wb, Ihc-P, Elisa - Q13535

Item Number: CM0013514



Introduction

Rabbit polyclonal antibody targeting ATR kinase involved in DNA damage checkpoint signaling. Suitable for WB, IHC-P, and ELISA; detects human, mouse, and rat ATR (~301 kDa). Ideal for cancer and genomic instability research.

[Learn More](#)

Parameter	Value
Product Name	ATR Rabbit pAb
Remarks/Alias	FRP1; MEC1; SCKL; FCTCS; SCKL1; ATR
Species	Human
GeneID (Human)	545
GeneID	545
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-500 of human ATR (NP_001175.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13535
Protein Weight	301kDa
Shipping	Ice bag