


Anti-Atr Rabbit Monoclonal Antibody For Wb And Elisa - Q13535

Item Number: CM0003527



Introduction

ATR Rabbit Monoclonal Antibody validated for Western Blot and ELISA, recognizing human and mouse ATR (301 kDa). Ideal for DNA damage response, checkpoint signaling, and nuclear envelope integrity research.

[Learn More](#)

Parameter	Value
Product Name	ATR Rabbit mAb
Remarks/Alias	FRP1; MEC1; SCKL; FCTCS; SCKL1; ATR
Species	Human
GeneID (Human)	545
GeneID	545
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 2100-2200 of human ATR (NP_001175.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13535
Protein Weight	301kDa
Shipping	Ice bag