


Anti-Atrip Monoclonal Antibody For Wb, Elisa - Q8Wxe1

Item Number: CM0003515



Introduction

Rabbit monoclonal antibody targeting human ATRIP (ATR-interacting protein), validated for WB and ELISA. Reacts with human, mouse, and rat samples. Ideal for DNA damage checkpoint signaling studies.

[Learn More](#)

Parameter	Value
Product Name	ATRIP Rabbit mAb
Remarks/Alias	ATRIP; ATR-interacting protein
Species	Human
Gene ID	84126
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human ATRIP (Q8WXE1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WXE1
Protein Weight	86kDa
Shipping	Ice bag