


Anti-Atrip Rabbit Polyclonal Antibody - Q8Wxe1

Item Number: CM0028411



Introduction

High-quality rabbit polyclonal antibody targeting human ATRIP (ATR-interacting protein), suitable for WB, IF/ICC, IP, and ELISA applications. Key for DNA damage response research.

[Learn More](#)

Parameter	Value
Product Name	ATRIP Rabbit pAb
Remarks/Alias	ATRIP
Species	Human
GeneID (Human)	84126
GeneID	84126
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human ATRIP (NP_569055.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human
SWISS	Q8WXE1
Protein Weight	86kDa
Shipping	Ice bag