

Anti-Atrnl1 Polyclonal Antibody For Wb And Elisa - Q5Vv63

Item Number: CM0016987

Introduction

High-quality anti-ATRNL1 rabbit polyclonal antibody for WB and ELISA, cross-reactive with human and mouse. Suitable for studying melanocortin signaling and energy homeostasis. Molecular weight 153kDa.

[Learn More](#)

Parameter	Value
Product Name	ATRNL1 Rabbit pAb
Remarks/Alias	ALP; bA338L11.1; bA454H24.1; ATRNL1
Species	Human
GeneID (Human)	26033
GeneID	26033
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 53-160 of human ATRNL1 (NP_001263211.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q5VV63
Protein Weight	153kDa
Shipping	Ice bag