


Anti-Atrx Rabbit Polyclonal Antibody For If-P And Elisa - P46100

Item Number: CM0022444



Introduction

Rabbit polyclonal antibody against human ATRX, suitable for immunofluorescence (IF-P) and ELISA. Cross-reacts with mouse. Ideal for chromatin remodeling and telomere biology research.

[Learn More](#)

Parameter	Value
Product Name	ATRX Rabbit pAb
Remarks/Alias	JMS; XH2; XNP; MRX52; RAD54; RAD54L; ZNF-HX
Species	Human
GeneID (Human)	546
GeneID	546
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2273-2413 of human ATRX (NP_000480.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P, ELISA
Cross Reactivity	Mouse
SWISS	P46100
Protein Weight	283kDa
Shipping	Ice bag