


Anti-Atrx Rabbit Polyclonal Antibody For Wb, Elisa - P46100

Item Number: CM0022010



Introduction

Rabbit polyclonal antibody against human ATRX, a chromatin remodeler involved in telomere maintenance and transcriptional regulation. Validated for WB and ELISA; cross-reacts with mouse ATRX. Suitable for epigenetic and oncology research.

[Learn More](#)

Parameter	Value
Product Name	ATRX Rabbit pAb
Remarks/Alias	JMS, XH2, XNP, MRX52, RAD54, RAD54L, ZNF-HX, ATRX
Species	Human
Gene ID (Human)	546
Gene ID	546
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 2100-2200 of human ATRX (NP_000480.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P46100
Protein Weight	283kDa
Shipping	Ice bag