

Anti-Atxn10 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ubb4

Item Number: CM0028287

Introduction

ATXN10 Rabbit Polyclonal Antibody for Western blot, immunofluorescence, and ELISA. Cross-reacts with human, mouse, rat. Ideal for studying spinocerebellar ataxia, cytokinesis, and neuronal survival.

[Learn More](#)

Parameter	Specification
Product Name	ATXN10 Rabbit pAb
Remarks/Alias	E46L; ATX10; SCA10; HUMEEP; ATXN10
Species	Human
Gene ID (Human)	25814
Gene ID	25814
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 346-475 of human ATXN10 (NP_037368.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBB4
Protein Weight	53kDa
Shipping	Ice bag