


Anti-Atxn2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q99700

Item Number: CM0022978



Introduction

Rabbit polyclonal antibody against human ATXN2 (Ataxin-2). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for neurodegeneration and EGFR trafficking research.

[Learn More](#)

Parameter	Value
Product Name	ATXN2 Rabbit pAb
Remarks/Alias	ATX2; SCA2; TNRC13; ATXN2
Species	Human
Gene ID (Human)	6311
Gene ID	6311
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1250-1331 of human ATXN2 (NP_002964.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99700
Protein Weight	140kDa
Shipping	Ice bag