


Anti-Auh Monoclonal Antibody For Wb And Elisa - Q13825

Item Number: CM0007050



Introduction

AUH Rabbit mAb (Q13825) validated for WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Suitable for studying mitochondrial leucine catabolism and RNA-binding functions.

[Learn More](#)

Parameter	Value
Product Name	AUH Rabbit mAb
Remarks/Alias	AUH
Species	Human
GeneID (Human)	549
GeneID	549
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 68-339 of human AUH(NP_001689.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13825
Protein Weight	36kDa
Shipping	Ice bag