

Anti-Auh Rabbit Polyclonal Antibody For Wb, Elisa - Q13825

Item Number: CM0026835



Introduction

High-quality anti-AUH rabbit polyclonal antibody validated for Western blot (WB) and ELISA. Reacts with human, mouse, and rat. Targets methylglutaconyl-CoA hydratase involved in leucine degradation and mitochondrial metabolism. Ideal for research in metabolic disorders.

[Learn More](#)

Parameter	Value
Product Name	AUH Rabbit pAb
Remarks/Alias	AUH
Species	Human
GeneID (Human)	549
GeneID	549
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 68-339 of human AUH (NP_001689.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13825
Protein Weight	36kDa
Shipping	Ice bag