


Anti-Aurka Polyclonal Antibody For Wb And Elisa - O14965

Item Number: CM0018025



Introduction

This rabbit polyclonal antibody targets human Aurora kinase A (AURKA), a key mitotic kinase. Validated for Western blot and ELISA, it detects a ~46 kDa band using a recombinant immunogen (1-130 aa). Suitable for oncology and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	AURKA Rabbit pAb
Remarks/Alias	AIK; ARK1; AURA; BTAK; STK6; STK7; STK15; PPP1R47; AURKA
Species	Human
GeneID (Human)	6790
Immunogen	Recombinant Protein / Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human AURKA (NP_940837.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14965
Protein Weight	46kDa
Shipping	Ice bag