


Anti-Aurkc Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Uqb9

Item Number: CM0026852



Introduction

Rabbit polyclonal antibody against human AURKC (Aurora kinase C), validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for studying chromosomal passenger complex functions, histone H3 phosphorylation, and cell division.

[Learn More](#)

Parameter	Value
Product Name	AURKC Rabbit pAb
Remarks/Alias	AIE2; AIK3; ARK3; AurC; SPGF5; STK13; HEL-S-90; aurora-C; AURKC
Species	Human
GeneID (Human)	6795
GeneID	6795
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human AURKC (NP_001015878.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UQB9
Protein Weight	36kDa
Shipping	Ice bag