

Anti-Aurora B Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96Gd4

Item Number: CM0008056



Introduction

Mouse monoclonal antibody against human Aurora B, applicable in WB, IHC-P, IF/ICC, and ELISA. Cross-reactive with mouse and rat. For cell division and mitosis research.

[Learn More](#)

Property	Value
Product Name	Aurora B Mouse mAb
Remarks/Alias	AIK2; AIM1; ARK2; AurB; IPL1; STK5; AIM-1; ARK-2; STK-1; STK12; PPP1R48; aurkb-sv1; aurkb-sv2
Species	Human
GeneID (human)	9212
GeneID	9212
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids1-100 of human TLR4 (NP_004208.2).
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96GD4
Protein Weight	39kDa
Shipping	Ice bag