

Anti-Auts2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Wxx7

Item Number: CM0020081

Introduction

Rabbit polyclonal antibody against AUTS2 (autism susceptibility gene 2 protein) validated for Western blot, IF/ICC, and ELISA. Detects human, mouse, and rat AUTS2. Suitable for chromatin remodeling and neurodevelopment studies.

[Learn More](#)

Parameter	Value
Product Name	AUTS2 Rabbit pAb
Remarks/Alias	MRD26, FBRSL2, AUTS2
Species	Human
Gene ID	26053
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 550-770 of human AUTS2 (NP_001120703).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WXX7
Protein Weight	139 kDa
Shipping	Ice bag