


Anti-Avil (Advillin) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O75366

Item Number: CM0019203



Introduction

Rabbit polyclonal antibody against human AVIL (Advillin), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying actin bundling, neuronal morphogenesis, ciliogenesis, and lamellipodia formation.

[Learn More](#)

Parameter	Value
Product Name	AVIL Rabbit pAb
Remarks/Alias	p92; DOC6; ADVIL; NPHS21; AVIL
Species	Human
GeneID (human)	10677
GeneID	10677
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 470-819 of human AVIL (NP_006567.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75366
Protein Weight	92kDa
Shipping	Ice bag