

Anti-Avp Polyclonal Antibody For Wb, Elisa - P01185

Item Number: CM0022927



Introduction

Rabbit polyclonal antibody targeting human AVP (vasopressin-neurophysin 2-copeptin) for WB and ELISA. Validated in human, mouse, and rat samples. Ideal for neuroendocrine research.

[Learn More](#)

Parameter	Value
Product Name	AVP Rabbit pAb
Remarks/Alias	VP; ADH; ARVP; AVRP; AVP-NP11; AVP
Species	Human
GeneID (Human)	551
GeneID	551
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01185
Protein Weight	17kDa
Shipping	Ice bag