

Anti-Avpi1 Polyclonal Antibody For Wb And Elisa - Q5T686

Item Number: CM0017924

Introduction

Rabbit polyclonal antibody against AVPI1 (VIP32/VIT32), validated for human, mouse, and rat samples in WB and ELISA. Suitable for MAP kinase, ENaC, and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	AVPI1 Rabbit pAb
Remarks/Alias	VIP32; VIT32; PP5395; AVPI1
Species	Human
Gene ID (Human)	60370
Gene ID	60370
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-147 of human AVPI1 (NP_068378.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5T686
Protein Weight	17kDa
Shipping	Ice bag