

Anti-Avpr1A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P37288

Item Number: CM0026715



Introduction

Rabbit polyclonal antibody against human AVPR1A (Vasopressin V1a receptor) for WB, IF/ICC, ELISA. Validated in human, mouse, rat. Recombinant immunogen (aa 349-418). Suitable for social behavior research.

[Learn More](#)

Parameter	Value
Product Name	AVPR1A Rabbit pAb
Remarks/Alias	V1aR; AVPR1; AVPR V1a; AVPR1A
Species	Human
Gene ID	552
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 349-418 of human AVPR1A (NP_000697.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37288
Protein Weight	47kDa
Shipping	Ice bag