

Anti-Avpr2 Rabbit Polyclonal Antibody For Wb, Elisa - P30518

Item Number: CM0026865



Introduction

Recombinant protein-based rabbit polyclonal antibody targeting human vasopressin V2 receptor (AVPR2). Validated for WB and ELISA. Ideal for renal water reabsorption and G protein-coupled receptor studies.

[Learn More](#)

Parameter	Value
Product Name	AVPR2 Rabbit pAb
Remarks/Alias	DI1; DIR; NDI; V2R; ADHR; DIR3; NDI1; AVPR2
Species	Human
GeneID (Human)	554
GeneID	554
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 232-371 of human AVPR2 (NP_000045.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P30518
Protein Weight	40kDa
Shipping	Ice bag