

Anti-Axin1 Rabbit Polyclonal Antibody For Wb And Elisa - O15169

Item Number: CM0014638



Introduction

High-quality anti-AXIN1 rabbit polyclonal antibody validated for Western blotting and ELISA, with cross-reactivity to human, mouse, and rat AXIN1. Ideal for studying Wnt/ β -catenin signaling and tumor suppressor mechanisms.

[Learn More](#)

Parameter	Value
Product Name	AXIN1 Rabbit pAb
Remarks/Alias	AXIN; PPP1R49; AXIN1
Species	Human
GeneID (Human)	8312
GeneID	8312
Immunogen	Recombinant fusion protein corresponding to amino acids 1-351 of human AXIN1 (NP_851393.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15169
Protein Weight	96kDa
Shipping	Ice bag