

Anti-Axin2 Polyclonal Antibody For Wb, Elisa - Q9Y2T1

Item Number: CM0019468



Introduction

Rabbit polyclonal antibody targeting human AXIN2 (Axin-2/Conductin) validated for western blot and ELISA; cross-reacts with mouse and rat. An essential tool for Wnt signaling pathway research.

[Learn More](#)

Parameter	Value
Product Name	AXIN2 Rabbit pAb
Remarks/Alias	AXIL; ODCRC5; AXIN2
Species	Human
Gene ID (human)	8313
Gene ID	8313
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 400-500 of human AXIN2 (NP_004646.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9Y2T1
Protein Weight	94kDa
Shipping	Ice bag