

Anti-Axin2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y2T1

Item Number: CM0020369



Introduction

Rabbit polyclonal antibody against human AXIN2, validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, rat. Ideal for Wnt signaling pathway studies. SWISS: Q9Y2T1.

[Learn More](#)

Parameter	Value
Product Name	AXIN2 Rabbit pAb
Remarks/Alias	AXIL; ODCRC5; AXIN2
Species	Human
Gene ID (Human)	8313
Gene ID	8313
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 743-843 of human AXIN2 (NP_004646.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2T1
Protein Weight	94kDa
Shipping	Ice bag