


Anti-Axl Monoclonal Antibody For Wb - P30530

Item Number: CM0004767



Introduction

Rabbit monoclonal antibody targeting human Axl receptor tyrosine kinase (UFO) for Western blot and ELISA applications. Reacts with human and mouse Axl, suitable for cell signaling and oncology studies.

[Learn More](#)

Parameter	Value
Product Name	Axl Rabbit mAb
Remarks/Alias	ARK; UFO; JTK11; Tyro7; Axl
Species	Human
GeneID (Human)	558
GeneID	558
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-449 of human Axl (NP_068713.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P30530
Protein Weight	98kDa
Shipping	Ice bag