


Anti-Axl Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P30530

Item Number: CM0021064



Introduction

Rabbit polyclonal antibody against human AXL (UFO) tyrosine kinase receptor, validated for WB, IHC-P, IF/ICC, IP, and ELISA. Suitable for cancer and viral entry research.

[Learn More](#)

Parameter	Value
Product Name	Axl Rabbit pAb
Remarks/Alias	ARK; UFO; JTK11; Tyro7; Axl
Species	Human
GeneID (Human)	558
GeneID	558
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-449 of human AXL (NP_068713.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human
SWISS	P30530
Protein Weight	98kDa
Shipping	Ice bag