


Anti-Azi2 Rabbit Monoclonal Antibody For Wb And Elisa - Q9H6S1

Item Number: CM0003474



Introduction

High-quality rabbit monoclonal AZI2 antibody validates for Western blotting and ELISA. Recognizes human, mouse, and rat AZI2 (NAP1/TILP), a key adaptor in antiviral innate immunity and NF- κ B signaling. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	AZI2 Rabbit mAb
Remarks/Alias	AZ2, NAP1, TILP, AZI2
Species	Human
GeneID (Human)	64343
GeneID	64343
Immunogen	Recombinant protein containing a sequence corresponding to amino acids 282-387 of human AZI2 (NP_071906.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H6S1
Protein Weight	45kDa
Shipping	Ice bag