

Anti-B-Raf Polyclonal Antibody For Wb, Elisa - P15056

Item Number: CM0030894



Introduction

Rabbit polyclonal antibody against human B-Raf, validated for WB and ELISA. Detects B-Raf in human, mouse, and rat samples. Ideal for MAP kinase signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	B-Raf Rabbit pAb
Remarks/Alias	NS7; B-raf; BRAF1; RAFB1; B-RAF1; BRAF-1; B-Raf
Species	Human
Gene ID (Human)	673
Gene ID	673
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 153-313 of human B-Raf (NP_004324.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15056
Protein Weight	84kDa
Shipping	Ice bag