

Anti-B-Raf Rabbit Pab For Wb, Ihc-P, Elisa - P15056

Item Number: CM0014747



Introduction

High-quality rabbit polyclonal antibody against human B-Raf, validated for WB, IHC-P, ELISA. Suitable for studying MAPK signaling in cancer research. Cross-reacts with mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	B-Raf Rabbit pAb
Remarks/Alias	NS7; B-raf; BRAF1; RAFB1; B-RAF1; BRAF-1; B-Raf
Species	Human
GeneID (human)	673
GeneID	673
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 350-449 of human B-Raf (NP_004324.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15056
Protein Weight	84kDa
Shipping	Ice bag