

Anti-B-Raf V600E Polyclonal Antibody For Wb And Elisa - P15056

Item Number: CM0020829



Introduction

High-specificity rabbit polyclonal antibody recognizing the oncogenic B-Raf V600E mutant. Validated for Western blotting and ELISA; ideal for MAPK pathway signaling studies and cancer research.

[Learn More](#)

Parameter	Value
Product Name	B-Raf V600E Rabbit pAb
Remarks/Alias	NS7; B-raf; BRAF1; RAFB1; B-RAF1; BRAF-1; B-Raf
Species	Human
GeneID (Human)	673
GeneID	673
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 507-606 of B-Raf V600E (NP_004324.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P15056
Protein Weight	84kDa
Shipping	Ice bag