

Anti-B3Galt4 Rabbit Polyclonal Antibody For Wb, Elisa - 096024

Item Number: CM0018403

Introduction

Rabbit polyclonal antibody targeting human B3GALT4, a Golgi-resident galactosyltransferase involved in ganglioside biosynthesis. Validated for WB and ELISA; cross-reacts with mouse. Ideal for diagnostic research and glycobiology studies.

[Learn More](#)

Parameter	Value
Product Name	B3GALT4 Rabbit pAb
Remarks/Alias	GALT2; GALT4; BETA3GALT4; B3GALT4
Species	Human
GeneID (Human)	8705
GeneID	8705
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-270 of human B3GALT4 (NP_003773.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O96024
Protein Weight	42kDa
Shipping	Ice bag