

Anti-B3Gnt6 Polyclonal Antibody For Wb, Elisa - Q6Zmb0

Item Number: CM0012996

Introduction

High-quality rabbit polyclonal antibody targeting human B3GNT6 (core 3 synthase). Suitable for Western blot and ELISA, with reported cross-reactivity to rat. Ideal for mucin-type O-glycosylation and gastrointestinal research.

[Learn More](#)

Parameter	Value
Product Name	B3GNT6 Rabbit pAb
Remarks/Alias	BGnT-6; B3Gn-T6; beta3Gn-T6; beta-1; 3-Gn-T6; B3GNT6
Species	Human
Gene ID (Human)	192134
Gene ID	192134
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-200 of human B3GNT6 (NP_619651.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q6ZMB0
Protein Weight	43kDa
Shipping	Ice bag