


Anti-B4Galt1 Rabbit Polyclonal Antibody For Wb And Elisa - P15291

Item Number: CM0021511



Introduction

Rabbit polyclonal antibody against human B4GALT1 for Western blot and ELISA. B4GALT1 is a Golgi-membrane galactosyltransferase involved in N-glycan synthesis and cell adhesion.

[Learn More](#)

Parameter	Specification
Product Name	B4GALT1 Rabbit pAb
Remarks/Alias	GT1; GTB; CDG2D; GGTB2; B4GAL-T1; beta4Gal-T1; B4GALT1
Species	Human
GeneID	2683 (Human)
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 50-215 of human B4GALT1 (NP_001488.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P15291
Protein Weight	44kDa
Shipping	Ice bag