


Anti-B4Galt1 Rabbit Polyclonal Antibody [Ko Validated] For Wb And Elisa - P15291

Item Number: CM0030485



Introduction

KO validated rabbit polyclonal antibody targeting human B4GALT1 (beta-1,4-galactosyltransferase 1). Suitable for Western blot (WB) and ELISA. Ideal for studies on glycosylation, cell adhesion, and lactose synthesis. UniProt P15291.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] B4GALT1 Rabbit pAb
Remarks/Alias	GT1; GTB; CDG2D; GGTB2; B4GAL-T1; beta4Gal-T1; T1
Species	Human
Gene ID (Human)	2683
Gene ID	2683
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 45-210 of human B4GALT1 (NP_001488.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P15291
Protein Weight	44kDa
Shipping	Ice bag