

## Anti-B4Galt2 Rabbit Polyclonal Antibody For Wb, Elisa - O60909

Item Number: CM0022645

### Introduction

Rabbit polyclonal antibody against human B4GALT2, validated for Western blot and ELISA. Reacts with human, mouse, and rat. Recombinant immunogen (aa 155-271). Predicted molecular weight 42 kDa.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | B4GALT2 Rabbit pAb                                                                                                    |
| Remarks/Alias    | B4Gal-T2; B4Gal-T3; beta4Gal-T2; B4GALT2                                                                              |
| Species          | Human                                                                                                                 |
| GeneID (Human)   | 8704                                                                                                                  |
| GeneID           | 8704                                                                                                                  |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 155-271 of human B4GALT2 (NP_085076.2). |
| Source           | Rabbit                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                 |
| Application      | WB, ELISA                                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                                     |
| SWISS            | O60909                                                                                                                |
| Protein Weight   | 42kDa                                                                                                                 |
| Shipping         | Ice bag                                                                                                               |