

Anti-B4Galt3 Polyclonal Antibody For Wb & Elisa - O60512

Item Number: CM0023118

Introduction

High-quality rabbit anti-B4GALT3 polyclonal antibody. Suitable for WB and ELISA applications. Detects endogenous levels of beta-1,4-galactosyltransferase 3 in human, mouse, and rat samples. Ideal for glycosylation and glycolipid metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	B4GALT3 Rabbit pAb
Remarks/Alias	beta4Gal-T3; B4GALT3
Species	Human
GeneID (Human)	8703
GeneID	8703
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-150 of human B4GALT3 (NP_003770.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60512
Protein Weight	44kDa
Shipping	Ice bag