

Anti-B4Galt4 Rabbit Polyclonal Antibody For Wb, Elisa - O60513

Item Number: CM0018191

Introduction

Rabbit polyclonal antibody against human B4GALT4 (beta-1,4-galactosyltransferase 4) for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for glycobiology and biomarker research.

[Learn More](#)

Parameter	Value
Product Name	B4GALT4 Rabbit pAb
Remarks/Alias	B4Gal-T4; beta4Gal-T4; B4GALT4
Species	Human
GeneID (Human)	8702
GeneID	8702
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-344 of human B4GALT4 (NP_003769.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60513
Protein Weight	40kDa
Shipping	Ice bag