

Anti-B9D1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Upm9

Item Number: CM0030021

Introduction

Rabbit polyclonal antibody against human B9D1, validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for cilia and ciliopathy research.

[Learn More](#)

Parameter	Specification
Product Name	B9D1 Rabbit pAb
Remarks/Alias	B9; MKS9; EPPB9; MKSR1; JBTS27; MKSR-1; B9D1
Species	Human
GeneID (Human)	27077
GeneID	27077
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human B9D1 (NP_056496.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UPM9
Protein Weight	23kDa
Shipping	Ice bag