

Anti-BaalC Rabbit Polyclonal Antibody For Wb, Elisa - Q8Wxs3

Item Number: CM0013685

Introduction

High-quality rabbit polyclonal antibody against BAALC, validated for WB and ELISA. Cross-reacts with human, mouse, and rat BAALC. Supports research in acute leukemia and neurobiology.

[Learn More](#)

Parameter	Value
Product Name	BAALC Rabbit pAb
Remarks/Alias	BAALC
Species	Human
GeneID (Human)	79870
GeneID	79870
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-145 of human BAALC (NP_079088.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WXS3
Protein Weight	16kDa
Shipping	Ice bag