

Anti-Baat Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q14032

Item Number: CM0029522

Introduction

Rabbit polyclonal antibody against human BAAT (bile acid-CoA:amino acid N-acyltransferase). Validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Immunogen: recombinant protein (aa 309-418).

[Learn More](#)

Parameter	Value
Product Name	BAAT Rabbit pAb
Remarks/Alias	BAT; HCHO; BACAT; BACD1; BAAT
Species	Human
GenelD (Human)	570
GenelD	570
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 309-418 of human BAAT (NP_001692.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14032
Protein Weight	46kDa
Shipping	Ice bag