


## Anti-Bace1 Rabbit Monoclonal Antibody For Wb And Elisa - P56817

Item Number: CM0003829



### Introduction

Rabbit monoclonal antibody targeting human BACE1 (beta-secretase 1), validated for WB and ELISA. Cross-reacts with mouse and rat. Essential for Alzheimer's disease research and diagnostics.

[Learn More](#)

| Parameter        | Value                                                                                                                   |
|------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Name     | BACE1 Rabbit mAb                                                                                                        |
| Remarks/Alias    | ASP2; BACE; HSPC104; BACE1                                                                                              |
| Species          | Human                                                                                                                   |
| GeneID (Human)   | 23621                                                                                                                   |
| GeneID           | 23621                                                                                                                   |
| Immunogen        | Synthetic Peptide   A synthetic peptide corresponding to a sequence within amino acids 402-501 of human BACE1 (P56817). |
| Source           | Rabbit                                                                                                                  |
| Category         | Monoclonal Antibodies                                                                                                   |
| Application      | WB, ELISA                                                                                                               |
| Cross Reactivity | Human, Mouse, Rat                                                                                                       |
| SWISS            | P56817                                                                                                                  |
| Protein Weight   | 56kDa                                                                                                                   |
| Shipping         | Ice bag                                                                                                                 |