


Anti-Bach1 Rabbit Monoclonal Antibody For Wb, Elisa, Chip - O14867

Item Number: CM0006345



Introduction

Recombinant rabbit monoclonal antibody targeting transcription regulator BACH1 (SWISS O14867), validated for WB, ELISA, CHIP. Cross-reacts with human and mouse. Ideal for investigating oxidative stress and gene regulation pathways.

[Learn More](#)

Parameter	Value
Product Name	BACH1 Rabbit mAb
Remarks/Alias	BACH-1; BTBD24; BACH1
Species	Human
GeneID (human)	571
GeneID	571
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 410-736 of human BACH1 (NP_001177.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA, CHIP
Cross Reactivity	Human, Mouse
SWISS	O14867
Protein Weight	82kDa
Shipping	Ice bag