


Anti-Bach1 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - O14867

Item Number: CM0025680



Introduction

Rabbit polyclonal antibody targeting BACH1 (BTB and CNC homolog 1), validated for WB, IP, ELISA; cross-reacts with human, mouse. Ideal for studying oxidative stress and NFE2L2 pathway regulation.

[Learn More](#)

Parameter	Value
Product Name	BACH1 Rabbit Polyclonal Antibody
Remarks/Alias	BACH-1; BTBD24; BACH1
Species	Human
Gene ID (Human)	571
Gene ID	571
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-320 of human BACH1 (NP_996749.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	O14867
Protein Weight	82kDa
Shipping	Ice bag