

Anti-Bad Monoclonal Antibody For Wb, Ihc-P, Elisa - Q92934

Item Number: CM0003415



Introduction

Rabbit monoclonal antibody against human Bad (Bcl2-associated agonist of cell death), suitable for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	Bad Rabbit mAb
Remarks/Alias	BBC2; BCL2L8; Bad
Species	Human
GeneID (Human)	572
GeneID	572
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Bad (Q92934))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92934
Protein Weight	18kDa
Shipping	Ice bag