


Anti-Bad Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q92934

Item Number: CM0014862



Introduction

Rabbit polyclonal antibody targeting human Bad protein. Validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for apoptosis studies and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Bad Rabbit pAb
Remarks/Alias	BBC2; BCL2L8; Bad
Species	Human
Gene ID (Human)	572
Gene ID	572
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human Bad (NP_004313.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92934
Protein Weight	18kDa
Shipping	Ice bag