

Anti-Baf57/Smarce1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q969G3

Item Number: CM0015404



Introduction

KO-validated rabbit polyclonal antibody against BAF57/SMARCE1, a SWI/SNF chromatin remodeling subunit. Recommended for WB, IF/ICC, IP, ELISA. Cross-reacts with human, mouse, rat. Immunogen: recombinant protein (aa 1-411).

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] BAF57/SMARCE1 Rabbit pAb
Remarks/Alias	CSS5; BAF57; E1
Species	Human
GeneID (Human)	6605
GeneID	6605
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-411 of human BAF57/BAF57/SMARCE1 (NP_003070.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q969G3
Protein Weight	47kDa
Shipping	Ice bag