


## Anti-Baf60A Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96Gm5

Item Number: CM0028323



### Introduction

Rabbit polyclonal antibody targeting BAF60a (SMARCD1), validated for WB, IHC-P, and ELISA. Reacts with human, mouse, rat. Ideal for chromatin remodeling and neurodevelopment research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | BAF60a Rabbit pAb                                                                                                  |
| Remarks/Alias    | Baf60a; D15Kz1; BAF60a                                                                                             |
| Species          | Mouse                                                                                                              |
| GeneID (Human)   | 6602                                                                                                               |
| GeneID           | 83797                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of mouse BAF60a (NP_114030.2). |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, ELISA                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q96GM5                                                                                                             |
| Protein Weight   | 58kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |